

PACHPEN

Equipped with



*Precision at Your
Fingertips*

Accurate - Achieve improved patient outcomes with Accutome's revolutionary digital technology.

IOP Correction Calculation - Converts IOP measurements in seconds.

Easy to Use - Superior ergonomic design fits comfortably in any hand.

Versatile - Helps to save time and ease your practice demands.

Accutome

263 Great Valley Parkway
Malvern, PA 19355 USA
800-979-2020 610-889-0200
fax 610-889-3233
www.accutome.com

The Difference is Visible.



THE NEXT GENERATION IN PACHYMETRY

Gentle Touch Probe
ensures patient
comfort

Stores and averages
9 measurements

**NO NEED TO
RECHARGE -
Long-lasting
Lithium Battery**
will last approx.
15,000 readings

Range of 300 to 999 microns
Accuracy of +/- 5 microns
Resolution of +/- 1 micron
65 MHz sampling probe
Lightweight - 3 oz. (85 g.)

Accutome

263 Great Valley Parkway
Malvern, PA 19355 USA
800-979-2020 610-889-0200
Fax 610-889-3233
www.accutome.com
eyeopener@accutome.com



-PachPen®
Pachymeter