


The Ultimate Ophthalmic Examination Unit

Ideal for Wheelchair access

- Electric rise and fall on table
- Electric rise and fall on chair
- Fixed back or reclinable back chair
- Stand alone or combination unit
- Table locks in any position via strong magnetic lock
- Made in Milan of steel and wood



This Italian manufactured Duo examination unit is as robust as they come. Completely made of steel and wood this sturdy attractive unit comes equipped with electric rise and fall table top.

The Duo unit also has a rotating table with a strong electromagnetic lock enabling the user to position the table top at any angle making it ideal for wheelchair bound patients.

The unit is available as a stand- alone unit (pictured) or combined chair and stand unit.

Various combinations of non-reclining or reclining back chairs are available.

Choose from extensive colour schemes or simply have us match your decor! The stand pictured features aluminium panelling with dark grey trim. Chair is reclinable and features flip up arm rests.

Table height and instruments are controlled by switches on the table top and panel, however even this is customisable.