

OcuLight® GL Laser

High Performance and Value

- Superb efficiency with a full 1200 mW of deliverable laser energy
- Provides pulse durations and pulse intervals from 50 to 1000 milliseconds
- Expands system performance with convenient accessories
 - Includes a full line of delivery devices
 - Offers a broad selection of EndoProbe® handpieces for a variety of surgical procedures
 - Offers an optional wireless footswitch to reduce clutter and improve ergonomics and control
- Offers excellent ergonomics and true portability
 - Quiet operation eliminates distraction during treatment
 - Convenient, easily understood controls permit smooth interface with the laser system and delivery devices
 - Transports easily between treatment rooms and offices



IRIDEX

Elegantly simple solutions™

OcuLight® GL Laser

Specifications

Wavelength:	532 nm Green
Weight:	16.4 lb (7.5 kg)
Dimensions:	6" H × 12" W × 12" D (15 cm × 30 cm × 30 cm)
Connector Type:	Resistor
Electrical:	100–240 VAC, 50/60 Hz
Cooling:	Whisper fan
Exposure Duration:	CW-Pulse™: 50–1000 ms
Exposure Interval:	CW-Pulse: 50–1000 ms
Aiming Laser:	Diode laser, 635 nm nominal
Delivery Device	Portable SLA: 0–1000 mW
Power Output:	LIO: 0–1000 mW EndoProbe®: 0–1200 mW

CE 0086

Specifications are subject to change without notice. OcuLight, EndoProbe, IRIDEX and the IRIDEX logo are registered trademarks and TruFocus and CW-Pulse are trademarks of IRIDEX Corporation.

Products are covered by one or more of the following U.S. patents: *5,372,595; 5,511,085; 5,521,932; 5,663,979; 5,982,789; 5,979,554; 6,141,143; 6,144,484; 6,222,869; 6,327,291; 6,377,599; 6,540,391; 6,733,490; 6,800,076; 7,537,593; 7,766,904; 7,771,417; and 7,909,816.*
Clinical references available upon request.

EC REP Emergo Europe
Molenstraat 15, 2513 BH The Hague, The Netherlands, Tel.: (31) (0) 70 345-8570, Fax: (31) (0) 70 346-7299



IRIDEX

Elegantly simple solutions™

