

Ocular Saxena Retinal Grid 428 & 520



Product Code

OI-SRG428
OI-SRG520



*Designed with:
Sandeep Saxena, M.D.
Lucknow, India*

Design

- Monofilament line at 4.00mm spacing for Ocular 28D & 5.20mm for Ocular 20D spacing provides reference to the size of the optic disc.
- OI-SRG428 designed to fit anterior side of Ocular 28D Indirect lenses*
- OI-SRG520 designed to fit anterior side of Ocular 20D Indirect lenses*
- Estimate the size of inflammatory/non inflammatory retinal lesions.
- Grid spacing aides in estimating the size of ocular tumors.
- Easily estimate the amount of disk edema.
- OI-SRG428 ideal for retinopathy of prematurity (ROP)

*lens design with diamond knurl pattern only

Cleaning & Disinfecting

Per Cleaning Method 2



2255 116th Ave NE, Bellevue, Washington 98004-3039 USA
T: 425-455-5200 or 800-888-6616 F: 425-462-6669
E: ocular@ocular-instruments.com I: www.ocular-instruments.com

© 2006 OCULAR INSTRUMENTS
9742B3000